



OPTION PREPARATION BOX - B (PB-B)

The Preparation Box B is designed for preparation of sterile consumables in a GMP class B environment. The Preparation Box B is an optional solution that is placed next to a Dispensing Hot Cell or Shielded Glove Box. A shielded airtight transfer door to the adjacent cell is part of the Preparation Box B.

Several options such as outlets, cable entries and a laptop support are available for the PB.

| Variant | нст |
|---|---|
| Outside dimensions (mm) (W*D*H) | ± 850*1.330*2.650 |
| Inside dimensions (mm) (W*D*H) | ± 709*719*787 |
| Frontdoor opening (W*H) | ± 590*715 mm / 24*29 inch |
| Shielded transfer hatch to adjacent hot cell, opening [W*H] | ± 200*250 mm |
| Air inlet | With HEPA filter, including valve |
| Air outlet | Including valve |
| Air exhaust Air outlet | Max 30 m³/h |
| Power connection in accordance with local requirements | |
| Exterior smooth finish | Traffic white (RAL 9016), easy to decontaminate |
| Total weight: | ± 530 kg |



Model PB-B









STANDARD FEATURES

Operator / product safety

- Exterior smooth finish easy to decontaminate.
- Internal stainless steel box 316L, optimal micro-surface treatment.
- Interior finish complies with pharmaceutical regulations for cleaning and validation.
- Air inlet: HEPA filter.
- Air refreshes 20 times per hour.
- GMP Grade B area.

Model PB-B

• Low voltage control panel for all controls.

Functions

- Hinged acrylic door with two gloves and inflatable seals to maintain airtight integrity.
- Analog gauge for displaying the internal pressure of the compartment.
- Compression block for custom Airtight entry of tubes and cables.
- LED internal light fixtures, switchable from control panel.
- Duplex electrical outlet inside the box, switched from control panel.
- Lead is not used for any structural support.
- HPV preparation.

OPTIONAL FEATURES

- Particle counter.
- Digital pressure gauge instead of analog gauge.
- Preparation box with two shielded side doors.
- Universal support (outlet, cable pass through, arm for laptop).
- Extended cover plates to specified height.
- Configurable space beneath main compartment:
 - Maintenance hatch.
 - HPLC / storage compartment.
 - Shielded generator vault.
 - Laminar downflow for GMP class A classification.





18-10-2024 10:08:44

